



Pressure Infusion Infu-Surg®

The latex free Infu-Surg® Pressure Infusion Bag from Ethox Medical is a safe, high-quality, reliable device for intravenous and intra-arterial infusion treatment, including A-line pressure monitoring. Designed for single patient use, the Infu-Surg may help prevent and control cross-contamination and hospital acquired infections. The Infu-Surg is a cost-effective, trusted product that medical professionals have relied on for decades.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Infu-Surg® Pressure Infusion Bag

- Pressure relief valve prevents over inflation (330 mmHg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300 mmHg)
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested
- Loads quickly and easily
- Latex free, single patient, non-sterile

Cat#	Size	Description	Box
4005	500cc	Infu-Surg	5
4005H	500cc	Infu-Surg with Hook	5
4010	1000cc	Infu-Surg	5
4010H	1000cc	Infu-Surg with Hook	5
4030	3000cc	Infu-Surg	5



Pressure Infusion Surgi-Press®

The latex free Surgi-Press® Pressure Infuser from Ethox Medical is an effective, durable device for proper and accurate infusion. Designed for long term use, the Surgi-Press features a spring driven piston to provide consistent and accurate infusion. It delivers continuous pressure without periodic pumping, thus eliminating the need for constant monitoring.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Surgi-Press® Pressure Infuser

- Easy to use; does not require special disposables, selective training or unique cleaning methods
- Solid-driven piston exerts continuous pressure on fluid bag with minimum variation of pressure
- No bulb to pump; no gauge to monitor
- No periodic pumping or adjustments needed
- Wide opening for rapid and simple fluid bag insertion
- Designed for effective, long-term use
- Solid case molded from glass fiber reinforced polymers
- Transparent cover for easy visualization
- Cleans and stores easily

Cat#	Size	Description	Box
4081	500cc	Surgi-Press	1
4085	500cc	Surgi-Press Hyperbaric Version	1