



# ER2000™


*Simply the most innovative and easily accessible emergency cart available today! Specifically designed to take seconds off your emergency response time.*

A true emergency cart unlike any other on the market. For your emergency department, standardize on the ER2000™ Series. State-of-the-art design that is attractive yet functional. Simple and easy to use, press the bright yellow "open" button located at the top center. The doors spring open allowing full access to the contents needed in any emergency situation. Available in RED or BLUE.



# “Designed for today’s Healthcare Professionals”

- Quick and simple to operate
- Clear visible drawers
- Easy access to life saving supplies



## Standard with each ER2000™ Cart

### Dimensions (includes casters, bumper and top)

Overall.....26.5"W x 32"D x 44.25"H  
 Open.....43"W x 32"D x 44.25"H

### Drawer trays

(1-4) .....14.7"W x 20.7"D x 2.5"H  
 (5).....14.7"W x 20.7"D x 5.5"H

### Door bins

(4 each).....7.5"W x 2.9"D x 2.8"H  
 (2 each).....7.2"W x 5.6"D x 2.6"H



ABH-1 Ambu Bag Hooks (2)



CBB-ER Cardiac Board and Bracket



ER-LBL Airway, Breathing, Circulation Labels



IV-3 IV Pole and Mounting



O-2 Oxygen Tank Brackets



PS-100 Plastic Breakaway Seals