

Phlebotomy/Specimen Collection Cart

FEATURES and BENEFITS

- These multipurpose carts are made out of a combination of steel, aluminum, and plastic to create a lightweight solution that is still very sturdy.
- Removable ABS plastic molded top.
- Includes gatelock for security.
- The 3 and 4 inch drawers act as tray holders for removable dividers. We include these customizable divider trays in each cart.
- Ergonomic push handles.
- 12 clear plastic tilt bins built into the front of the cart that conveniently show their contents for fast and easy access.
- Toss away waste into the pass-through sharps container that locks inside the cart or into the built-in tilt-out waste bin (see below).
- Also includes a glove box holder and side mounted paper holder.
- 5" Premium easy roll casters. (1) full swivel, (2) swivel w/ brake, (1) swivel w/ tracking.
- Stabilizer Base with high impact, non-marking plastic bumper to add stability and help keep cart from tipping.



MTWA-34696-WHT

Dimensions:

Overall: 27"W x 24"D x 45"H

Drawers:

- (1) 10.60" W x 16.5"D x 3"H
- (2) 10.60" W x 16.5"D x 4"H
- (3) 10.60" W x 16.5"D x 6"H
- (4) 10.60" W x 16.5"D x 9"H
- (5) 22"W x 16.5"D x 6"H

MSWA-3469-WHT

Dimensions:

Overall: 27"W x 24"D x 37.25"H

Drawers:

- (1) 10.60" W x 16.5"D x 3"H
- (2) 10.60" W x 16.5"D x 4"H
- (3) 10.60" W x 16.5"D x 6"H
- (4) 10.60" W x 16.5"D x 9"H




Unique lockable area to hold sharps container.



Built in drawer dividers for easy removal and cleaning.



Tilt out waste bin.



Paper holder & glove box holder.