

# ViroVac™

## Smoke Plume Evacuation System

### Features:

Compact design allows for easy placement on a shelf or counter and accommodates stack-ability-- perfect for limited space locations

Ultra quiet operation will assist in providing a safe operating environment without disruption

35 hour maximum filter life delivers unsurpassed value and reduces per procedure costs

3-port filter design puts the power of hose and accessory options in your hand and eliminates adapters

Filter ports are covered for additional safety

Simplicity reduces training needs and enhances utilization



*Combining simplicity and ultra-quiet operation in a compact package, the ViroVac™ Surgical Smoke Evacuation System exceeds expectations for value and efficacy. The unique 3-port filter offers a variety of options when choosing an accessory or tube set to meet the clinical need.*

**BUFFALO Filter**  
World Leader in Surgical Smoke Evacuation™

ViroVac Lit Rev. B.pdf

# ViroVac™

Smoke Plume Evacuation System

ViroVac™ is designed to be simple to use, flexible in its options for tubing and accessories, and sets new precedents for quiet operation. Providing effective surgical smoke evacuation that encompasses 4 layers of filtration and up to 35 hours of filter life, the ViroVac will deliver the performance necessary to effectively manage surgical smoke plume with unsurpassed value. Operate with a footpedal, control panel, or optional Remote Switch Activator.



Versatile three-port filter will accommodate 1/4" (6mm), 3/8" (9.5mm) and 7/8" (22mm) tubing diameters can be used in various procedures.



Hidden male and female ports allow for ease of tube attachment without risk of user contact.



The optional EZ Link Automatic Activation Device turns the Smoke Evacuator on and off with ESU activation, (part: EZLINK01).

## ViroVac™ Smoke Plume Evacuation System

(Filters and tubing sold separately).

100/120 Volt System  
 PART # VV120  
 120VAC ± 10% 50Hz/60Hz  
 Dimensions:  
 15.2 cm H x 27.9 cm W x 39.4 cm D  
 (6" H x 11" W x 15.5"D)  
 Weight: 4.5kgs/10lbs

220/240 Volt System  
 PART # VV220  
 230VAC ± 10% 50Hz/60Hz  
 Dimensions:  
 15.2 cm H x 27.9 cm W x 39.4 cm D  
 (6" H x 11" W x 15.5"D)  
 Weight: 4.5kgs/10lbs

ViroVac™ ViroSafe® Filter  
 PART #VS353  
 Maximum Filter Life: 35 hours  
 Dimensions:  
 24 cm H x 14 cm W x 11.5 cm D  
 (9.5"H x 5.5"W x 4.5"D)  
 Weight: 1.13kgs/2.5lbs

\*Dependent on flow setting



BUFFALO FILTER™ — A12634



World leader in surgical smoke evacuation™

Made in the USA • CE  
 © 2012 Buffalo Filter

### Buffalo Filter®

5900 Genesee Street  
 Lancaster, NY 14086  
 U.S.A.

Phone: 800.343.2324 (US)

Phone: +1.716.835.7000

Fax: 716.835.3414

Email: [smoke-info@buffalofilter.com](mailto:smoke-info@buffalofilter.com)

Web: [www.buffalofilter.com](http://www.buffalofilter.com)