

Buffalo Filter's®

visiClear™

Surgical Smoke Plume Evacuators

Features:

SafePort Technology™ – Buffalo Filter's unique patent pending sliding port closure design ensures that cross contamination is reduced.

Intuitive Color Touch Screen – Increases productivity and accuracy of medical personnel while requiring minimal training.

Whisper Technology™ – Our technologically advanced control features allow for ultra-quiet operation with unprecedented sound quality.

Sleek Shelf Design – Lies flat on a surgical boom shelf with isolation feet that do not require hard attachment. Designed to accommodate stacking.

State-of-the-Art Electronics –

Features a 32-bit microprocessor with embedded color LCD driver to optimize speed and power.



Combining Whisper Technology™ for ultra-quiet operation, a sleek shelf design and touch screen convenience, the VisiClear™ System delivers cutting-edge protection against the hazards of surgical plume in an acute care environment. Buffalo Filter's SafePort Technology™ closed-port design delivers unsurpassed safety. The unique three-port filter offers a variety of options when choosing an accessory or hose set to meet the clinical need.

Buffalo Filter's® VisiClear™ Surgical Smoke Plume Evacuator

VisiClear™ is designed to be our most advanced acute care offering for surgical smoke plume evacuation. With a low procedure cost, VisiClear™ was designed using cutting-edge technology to effectively capture and filter surgical plume and remove odor, particulates and other potentially hazardous by-products in the smokiest of procedures. The unique three-port, four-stage high performance ULPA filter features Buffalo Filter's SafePort Technology™. The closed port design safely enables removal of hose sets and ensures that the possibility of cross contamination is reduced to a minimum.



The illuminated touch screen offers a clear and interactive display, giving your staff unmatched control in any clinical configuration.



Buffalo Filter's SafePort Technology™ not only makes attaching hose sets and accessories easy, it also helps to prevent cross contamination during procedures.



The optional EZ Link™ Automatic Activation Device, (sold separately), automatically turns the VisiClear™ System on and off with ESU activation, greatly extending filter life.

Noise Levels:

@ minimum power < 43dBA
@ maximum power < 55dBA

Dimensions (H x W x L):

18.14cm x 35.4cm x 40.87cm
(7.14in x 13.94in 16.09in)

Flow Rate:

22mm (7/8") 30 CFM / (839LPM)

Weight:

Without filter installed: 14.5lbs (6.58kg)
With filter installed: 17.5 lbs (7.94 kg)

Static Suction:

150mmHg vacuum

Certifications:

UL, cUI, CE

Made in USA



Product Ordering Information:

Catalog Number:	
100/120 Volt System	VC120
20/240 Volt System	VC220
Quantity:	1 per carton

Filter Ordering Information

Catalog Number:	VS135
-----------------	-------

Buffalo Filter®

5900 Genesee Street
Lancaster, NY 14086
U.S.A.

Phone: 800.343.2324 (US)

Phone: +1.716.835.7000

Fax: 716.835.3414

Email: smoke-info@buffalofilter.com

Web: www.buffalofilter.com

BF P/N VisiClear RevA

